



‘인알못’을 위한 심전도 인공지능 소개

Introduction of ECG AI for Beginners

11th Dec 2021

Mineok Chang, MD, PhD
VUNO Inc. Medical Director



Registered in R.O.K. & U.S. Patent and Trademark Office

Copyright Notice

The contents of these documents, including (but not limited to) all written material, images, and photos are copyright of VUNO, Inc. Unauthorized copying, reproduction, modification, republishing, transmission or distribution of any material without the prior written content of VUNO, Inc. constitute an infringement of intellectual property rights and may result in a criminal or civil action in accordance with Article 98 of Copyright Act. Copyright (C) VUNO Inc. All rights reserved.

Contents

What is AI ?


Qch: P 34
QRS 12
T 26
1.29

KR 557 MC
P 88 MC
PR 128 MC
QRS 72 MC
QT 314 MC
QTc 424 MC


P (II) 0.86 MB
S (V1) -0.79 MB
R (V5) 0.98 MB
Sokol. 1.87 MB



AI ? ML ? DL ?





Types of Machine Learning



Deep Learning

Artificial Neural Network (ANN) \cong Neural Network (NN)







By Geoffrey Hinton

Designing a Learning System = Modeling


$$y = ax + b$$

“ Training – Validation (Development) – Test ”

1. Problem Description
2. Choosing the Data (Parameter)
3. Choosing the Target Function
4. Choosing the Learning Architecture
5. Choosing the Learning Algorithm
6. Final Design (Hyperparameter)




$$D_G^*(\mathbf{x}) = \frac{p_{data}(\mathbf{x})}{p_{data}(\mathbf{x}) + p_g(\mathbf{x})}$$




ML in In Medicine

“ Training – Validation (Development) – Test “



Single Center


***Internal
Validation***



Multi Center

***External
Validation***

Medical Image on PACS



Contents

ECG data & analysis

Qch:
P 34
QRS 12
T 26
1,29



KR 557 MC
P 88 MC
PR 128 MC
QRS 72 MC
QT 314 MC
QTc 424 MC

P (II) 0.86 MB
S (V1) -0.79 MB
R (V5) 0.98 MB
Sokol. 1.87 MB



ECG is already Digitalized

What you see is PDF or JPEG on PACS or EMR
MUSE stores ECG raw data (XML, HL7)



Inference

x

XML data

```
CQAJAAKACAAIAAgACAAIAAgACAAIAAgACAAHAACACAAIAAgACAAIAAYABgAHAAGACAAH  
AAcACA[REDACTED]AgACAAIAAgABgAGAAcABWAGAAyABgAHAACABgAGAAyABgAGAAQA  
BAAFAAUABAAEAAQABQAFAAUABQFAAUABQAEAAQABAADAAMAawADAAIAAgACAAMBAAE  
AAQBAEAAEAAQAEAAQAAwACAAIAAgACAAIAAgACAAIAAwADAAMAawADAAIAAgADAAMA  
AgACAAIAAQABAAMAawADAAIAAgACAAIAAgACAAIAAgA[REDACTED]AEAAQADAAMAQAB  
AAEAAgACAAIAAgACAAIBAAEAAEA//8BAAMAawADAAMAawACAAIAAgACAAIAAgABAAEA  
BQAFAMAQAABAAMBAEAAEAAQCAAIABAAEAAAMAawADAAMAawADAAMAAGCAAQABAAD  
AAEAAQCAIAAgACAAIAAgACAAIAAgACAAMAawADAAIAAgACAAMAawACAAIAAgABAAEA  
AgACAAIAAgABAAAAAAAAAAAAAAAABAAMAQABAAEAAQAAAAAAAAAGCAAIAAgAGAAoEAAV  
ACEALAA7AE4YABxIEAiWCPAJQAmwCWAHsAXgBEADMAHQAFAPn/9//3//f/+//8CAAUA  
BwAHAAcABAD//wEAAQABAAAA/////wAAAAABAAEAAQABAAEAAQABAP//AQABAAIAAgAC  
AAEAAQABAAEAAgACAAEAAQABAAAAAAAABAIAAgACAAEAAQABAAIAAgADAAMAawADAAMA  
AwADAAIAAgACAAIAAgACAAIAAQABAAEABAAEAAIAAgACAAAAAABAAIAAgACAAIAAgAC  
AAMAawADAAMAawFAAUABQAEAAIAAgADAQAQABQFAAQABAAEAAQ[REDACTED]AFAAUA  
BQAFAAUABgAGAAyABQAFAAUACAAGAAQABgAGAAyABgAHAACACQAHAAKACwAIAAgACgAK  
AAoACgANAAsACwALAAoACgAMAawADAANAA4ADgAOAA4AEANAA0ADgAQABAADwAPABAA  
EAAQABAAEAAQABAA[REDACTED]ARABAAEAAQABAAEAAQAA8ADwARABIAEgARABEAEQAS  
ABIAEAAQABAAEAAQABAAEAAQABAAEAAQAA4ADgANA0ADAALAAsACwALAAsACwALAAKA  
CQAJAAKACQALAAwADAANAA4ADgAOAA4ACwALAAsACwAMAawACwALAAsACwALAAAsACgAI  
AAgACQAJAAoACgAJAAcACgAKAAoACQAIAGACAA[REDACTED]AgACQAJAAYABgAHAAKA  
CQAIAGACAAIAAgABgAGAAgACAAIAAgACAAIAAgACAAIAAYABgAIAAgACAAIAAgACQAJ  
AAgACAAIAAgACAAIAAgACAAIAAKACQAHAAcABwAIAAgACAAIAAgACAAIAAgACAAIAAgA  
CAAGAAyACAAIAAYABgAIAAgABwAHAACABwAIAAgACAAIAAgACAAIAAgACAAIAAgACgAK  
AAoACgALAAAsACwALAAAsACwALAAAsADAAMAawA
```

[</WaveFormData>](#)

$$y = ax + b$$

Model (Code)

```
def model_fn (params):  
  
    inputs = tf.keras.Input(name='waveform')  
  
    x = layers.Conv2D(  
        filters=64,  
        kernel_size=3,  
        strides=1,  
        padding='SAME',  
        use_bias=False,  
        kernel_initializer='he_normal')(x)  
    x = layers.BatchNormalization()(x)  
    x = layers.Activation('relu')(x)  
  
    x = block(64)(x)  
    x = block(64)(x)  
  
    x = block(128, upscale=True)(x)  
    x = block(128)(x)  
    x = block(128)(x)  
  
    x = block(256, upscale=True)(x)  
    x = block(256)(x)  
    x = block(256)(x)  
  
    x = layers.GlobalAveragePooling2D()(x)  
  
    model = tf.keras.Model(inputs=inputs, outputs=outputs)  
    model.summary()  
  
    model_spec = {}  
    model_spec['model'] = model  
    model_spec['loss'] = get_loss(params)  
    model_spec['optimizer'] = get_optimizer(params)  
  
    return x, model_spec
```

y

Output

$0 \sim 1$
Probability
With cut-off

MODEL_OUTPUT
0.655296028
0.542461872
0.644664586
0.902734876
0.801081657
0.576097906
0.660251319

Performance of Model

Confusion Matrix

	Actual Positive	Actual Negative
Model Positive	TP	FP (type I error)
Model Negative	FN (type 2 error)	TN

Sensitivity, Recall (재현율)	$\frac{TP}{TP + FN}$
Specificity	$\frac{TN}{FP + TN}$
PPV, Precision (정밀도)	$\frac{TP}{TP + FP}$
NPV	$\frac{TN}{TN + FN}$
F1 (조화평균)	$2 * \frac{Precision * Recall}{Precision + Recall} = \frac{2TP}{2TP + FP + FN}$
Accuracy	$\frac{TP + TN}{TP + FP + FN + TN}$

Data Imbalance matters

$$\text{Accuracy} = \frac{TP+TN}{TP+FP+FN+TN}$$

$$\text{F1 score} = \frac{2TP}{2TP+FP+FN}$$

Ex. Glioblastoma prevalence : 0.003%




Confusion Matrix		
	Actual Positive	Actual Negative
Model Positive	0 (TP)	0(FP)
Model Negative	3(FN)	100,000(TN)

$$\text{Glioblastoma prediction Accuracy} = \frac{100,000}{100,003} = 99.99\%$$


$$\text{Glioblastoma F1 score} = \frac{0}{3} = 0\%$$

Recall (Sensitivity) vs Precision (PPV)

Trade-off exist !
ROC is not sufficient !



Trade-off between Accuracy & Interpretability





"AI is Black Box"



Explainable AI, is it possible?

What you see is, what you know !



Contents

How ECG AI works?

Qch:
P 34
QRS 12
T 26
1,29

KR 557 MC
P 88 MC
PR 128 MC
QRS 72 MC
QT 314 MC
QTc 424 MC

P (II) 0.86 MB
S (V1) -0.79 MB
R (V5) 0.98 MB
Sokol. 1.87 MB




Amount of Information

Computer

vs.

Human




Amount of information is $12L \times 10\text{sec} = 120$
Considering augmented lead, total information would be **90**


Amount of information is **30**

More Data = For a feature → More accurate mean value


Normal Distribution



표본수:20개



표본수:200개




표본수:2,000개

Law of Large Numbers


표본집단의 크기가 커지면, 그 표본평균이 모평균에 가까워 짐. 즉, 취합하는 표본의 수가 많을 수록, 통계적 정확도는 올라간다.

More Data = More features


Rare feature could be retained.





표본수:20개



표본수:200개



표본수:2,000개



Empirical Risk Minimization : Goal of ML

$$R_{emp}(h) = \frac{1}{n} \sum_{i=1}^n L(h(x_i), Y_i)$$


XML data

Ground Truth


모델함수

평균


Loss Function : 모델과 Ground Truth와의 차이



표본수:20개




표본수:200개



표본수:2,000개

The more information, the better performance ?





Thank you

mineok.chang@vuno.co

oklizard81@gmail.com

T 02.515.6646

E hello@vuno.co

H www.vuno.co

VUNO® VUNOMed®

Registered in R.O.K. & U.S. Patent and Trademark Office

Copyright Notice

The contents of these documents, including (but not limited to) all written material, images, and photos are copyright of VUNO, Inc. Unauthorized copying, reproduction, modification, republishing, transmission or distribution of any material without the prior written content of VUNO, Inc. constitute an infringement of intellectual property rights and may result in a criminal or civil action in accordance with Article 98 of Copyright Act. Copyright (C) VUNO Inc. All rights reserved.